

CAS # 5614 – 37 - CYCLOPENTYL METHYL ETHER

Scientific OEM

208 Sujata Chambers , 1/3 Abhaychand Gandhi Marg, Masjid Bunder (W) , Mumbai - 400 009. India.

Tel : +91-22- 2343 7546 / 2341 3094

Fax : +91-22-2343 6572

E-MAIL: info@scientificoem.com WEB: www.scientificoem.com

1 - Product Information

Product Name CYCLOPENTYL METHYL ETHER

2 - Hazards Identification

SPECIAL INDICATION OF HAZARDS TO HUMANS AND THE ENVIRONMENT

Highly flammable. Harmful if swallowed. Irritating to eyes and skin.

3 - Composition/Information on Ingredients

Product Name	CAS #	EC no Annex I
Index Number		
CYCLOPENTYL METHYL ETHER	5614-37-9	None
None		

Formula C₆H₁₂O

Molecular Weight 100.16 AMU

4 - First Aid Measures

AFTER INHALATION

If inhaled, remove to fresh air. If not breathing give artificial respiration. If breathing is

difficult, give oxygen.

AFTER SKIN CONTACT

In case of contact, immediately wash skin with soap and copious amounts of water.

AFTER EYE CONTACT

In case of contact, immediately flush eyes with copious amounts of water for at least 15 minutes.

AFTER INGESTION

If swallowed, wash out mouth with water provided person is conscious. Call a physician.

5 - Fire Fighting Measures

EXTINGUISHING MEDIA

Suitable: For small (incipient) fires, use media such as "alcohol" foam, dry chemical, or carbon dioxide. For large fires, apply water from as far as possible. Use very large quantities (flooding) of water applied as a mist or spray; solid streams of water may be ineffective. Cool all affected containers with flooding quantities of water.

SPECIAL RISKS

Specific Hazard(s): Flammable liquid. Emits toxic fumes under fire conditions.

Explosion Hazards: Vapor may travel considerable distance to source of ignition and flash back. Container explosion may occur under fire conditions.

SPECIAL PROTECTIVE EQUIPMENT FOR FIREFIGHTERS

Wear self-contained breathing apparatus and protective clothing to prevent contact with skin and eyes.

6 - Accidental Release Measures

PERSONAL PRECAUTION PROCEDURES TO BE FOLLOWED IN CASE OF LEAK OR SPILL

Evacuate area. Shut off all sources of ignition. Use non sparking tools.

PROCEDURE(S) OF PERSONAL PRECAUTION(S)

Wear self-contained breathing apparatus, rubber boots, and heavy rubber gloves.

METHODS FOR CLEANING UP

Cover with dry-lime, sand, or soda ash. Place in covered containers using non-sparking tools and transport outdoors.

Ventilate area and wash spill site after material pickup is complete.

7 - Handling and Storage

HANDLING

Directions for Safe Handling: Do not breathe vapor. Do not get in eyes, on skin, on clothing. Avoid prolonged or repeated exposure.

STORAGE

Conditions of Storage: Keep tightly closed. Keep away from heat, sparks, and open flame.

SPECIAL REQUIREMENTS: Air and light sensitive. Handle and store under inert gas.

8 - Exposure Controls / Personal Protection

ENGINEERING CONTROLS

Safety shower and eye bath. Use nonsparking tools. Mechanical exhaust required.

GENERAL HYGIENE MEASURES

Remove and wash contaminated clothing promptly. Wash thoroughly after handling.

PERSONAL PROTECTIVE EQUIPMENT

Respiratory Protection: Use respirators and components tested and approved under appropriate government standards such as NIOSH (US) or CEN (EU). Where risk assessment shows air-purifying respirators are appropriate use a full-face respirator with multi-purpose combination (US) or type ABEK (EN 14387) respirator cartridges as a backup to engineering controls. If the respirator is the sole means of protection, use a full-face supplied air respirator.

Hand Protection: Compatible chemical-resistant gloves.

Eye Protection: Chemical safety goggles.

9 - Physical and Chemical Properties

Appearance Physical State: Clear liquid

Color: Colorless

Property Value At Temperature or Pressure

pH N/A

BP/BP Range 106 °C

MP/MP Range - 140.0 °C

Flash Point - 1.0 °C Method: closed cup

Flammability N/A

Autoignition Temp N/A

Oxidizing Properties N/A

Explosive Properties N/A

Explosion Limits Lower: 1.1 %

Upper: 9.9 %

Vapor Pressure N/A

SG/Density 0.86 g/cm³ 20 °C

Partition Coefficient Log Kow: 1.59 20 °C

Viscosity N/A

Vapor Density 3.5 g/l

Saturated Vapor Conc. N/A

Evaporation Rate N/A

Bulk Density N/A

Decomposition Temp. N/A

Solvent Content N/A

Water Content N/A

Surface Tension N/A

Conductivity N/A

Miscellaneous Data N/A

Solubility N/A

10 - Stability and Reactivity

STABILITY

Conditions to Avoid: Direct sunlight Potentially explosive peroxides can form on long time storage in contact with air. Light and heat accelerate peroxide formation. May form peroxides on contact with air.

Materials to Avoid: Oxidizing agents, Acids.

HAZARDOUS DECOMPOSITION PRODUCTS

Hazardous Decomposition Products: Carbon monoxide, Carbon dioxide.

HAZARDOUS POLYMERIZATION

Hazardous Polymerization: Will not occur

11 - Toxicological Information

ROUTE OF EXPOSURE

Skin Contact: Causes skin irritation.

Skin Absorption: May be harmful if absorbed through the skin.

Eye Contact: Causes eye irritation.

Inhalation: Material may be irritating to mucous membranes and upper respiratory tract.

May be harmful if inhaled.

Ingestion: Harmful if swallowed.

12 - Ecological Information

No data available.

ECOTOXICOLOGICAL EFFECTS

Test Type: LC0 Fish

Time: 96 h

Value: > 220 mg/l

13 - Disposal Considerations

SUBSTANCE DISPOSAL

Contact a licensed professional waste disposal service to dispose of this material. Burn in a chemical incinerator equipped with an afterburner and scrubber but exert extra care in igniting as this material is highly flammable. Observe all federal, state, and local environmental regulations.

14 - Transport Information

RID/ADR

UN#: 1993

Class: 3

PG: 2

Proper Shipping Name: Flammable liquid, n.o.s.

IMDG

UN#: 1993

Class: 3

PG: 2

Proper Shipping Name: Flammable liquid, n.o.s.

Marine Pollutant: No

Severe Marine Pollutant: No

Technical Name: Required

IATA

UN#: 1993

Class: 3

PG: 2

Proper Shipping Name: Flammable liquid, n.o.s.

Inhalation Packing Group I: No

Technical Name: Required

15 - Regulatory Information

CLASSIFICATION AND LABELING ACCORDING TO EU DIRECTIVES

INDICATION OF DANGER: F-Xn

Highly Flammable. Harmful.

R-PHRASES: 11-22-36/38

Highly flammable. Harmful if swallowed. Irritating to eyes and skin.

S-PHRASES: 26

In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.

Caution: Substance not yet fully tested (EU).

16 - Other Information

Prepared by: Scientific OEM

For questions contact: Scientific OEM

208 Sujata Chambers , 1/3 Abhaychand Gandhi Marg, Masjid Bunder (W) , Mumbai - 400 009. India.

Tel : +91-22- 2343 7546 / 2341 3094 Fax : +91-22-2343 6572 E-MAIL:

info@scientificoem.com

The statements contained herein are offered for informational purposes only and are based upon technical data that LOBA CHEMIE believes to be accurate. It is intended for use only by persons having the necessary technical skills and at their own discretion and risk. Since conditions and manner of use are outside our control, we make NO WARRANTY, EXPRESS OR IMPLIED, OF MERCHANTABILITY, and FITNESS OR OTHERWISE